

Product Information Sheet for NR-52948

SARS-Related Coronavirus 2, Wuhan-Hu-1 Spike-Pseudotyped Lentiviral Kit

Catalog No. NR-52948

For research use only. Not for human use.

Contributor:

Jesse Bloom, Associate Professor, Department of Basic Sciences, Fred Hutchinson Research Center, Seattle, Washington, USA and Alejandro Balasz, Assistant Professor, The Ragon Institute of MGH, MIT and Harvard, Cambridge, Massachusetts, USA

Manufacturer:

BEI Resources

Product Description:

NR-52948 is intended for producing pseudotyped particles/pseudovirions and is not for soluble protein expression.

The severe acute respiratory syndrome-related coronavirus 2 (SARS-CoV-2), Wuhan-Hu-1 (GenBank: [NC_045512](#)) spike-pseudotyped lentiviral kit (NR-52948) is designed to generate pseudotyped lentiviral particles with the spike (S) glycoprotein gene, as well as luciferase (Luc2) and green fluorescent protein (GFP). Protocols for the use of these items are published.^{1,2}

NR-52948 consists of the five plasmids listed in Table 1. Descriptions of each component are included below.

Table 1: SARS-CoV-2 Lentiviral Kit

Plasmid Type	Insert	BEI Resources Catalog Number
Viral Entry Protein	Spike Glycoprotein	NR-52514
Lentiviral Backbone	Luc2; ZsGreen	NR-52516
Helper Plasmid	Gag; pol	NR-52517
Helper Plasmid	Tat1b	NR-52518
Helper Plasmid	Rev1b	NR-52519

NR-52514 was designed by codon optimizing the S sequence for mammalian expression and subcloned into the pHDM vector under the cytomegalovirus (CMV) promoter.^{1,2}

NR-52516 was designed by fusing the synthetic firefly luciferase (Luc2) gene to the encephalomyocarditis internal ribosomal entry site (IRES) and synthetic *Zoanthus* sp. green fluorescent protein (ZsGreen1) gene, allowing simultaneous expression of Luc2 and the ZsGreen1 gene, which were subcloned into the pHAGE lentiviral backbone vector under the CMV promoter.^{1,2,3,4} The Luc2 gene has been codon optimized for mammalian expression and has had cryptic transcription factor binding sites removed. The ZsGreen1

gene has been codon optimized for mammalian expression and engineered for brighter fluorescence.³ In addition, the pHAGE vector includes the Woodchuck hepatitis virus post-transcriptional regulatory element to enhance levels of transcription and gene expression. The resulting size of the plasmid is approximately 9370 base pairs.

NR-52517 and NR-52518 were designed by codon optimizing the genes *gag* and *pol* and the gene *tat1b*, respectively, from the human immunodeficiency virus (HIV) and subcloning them into the pHDM vector under the CMV promoter.^{1,2,4} The resulting plasmid sizes are approximately 8910 base pairs and 4830 base pairs, respectively.

NR-52519 was designed by codon optimizing the *rev1b* gene from HIV and subcloning into the [pRC-CMV](#) vector under the CMV promoter.^{1,2,4} NR-52519 contains a neomycin (G418) selectable marker for mammalian expression. The resulting size of the plasmid is approximately 5900 base pairs.

NR-52514, NR-52516, NR-52517, NR-52518 and NR-52519 contain the beta-lactamase gene, *bla*, to provide transformant selection through ampicillin resistance in *Escherichia coli* (*E. coli*).

The complete plasmid sequences and maps are provided on the BEI Resources webpage. The plasmids were produced in *E. coli* and extracted.

Material Provided:

Each kit contains one vial of each plasmid DNA in 10 mM Tris-HCl, 1 mM EDTA, pH 8.0. The DNA concentrations and volumes provided are shown on the Certificate of Analysis. The vials should be centrifuged prior to opening. **Note:** The contents of each vial should be used to replicate the plasmid in *E. coli* prior to mammalian expression.

Packaging/Storage:

NR-52948 was packaged aseptically in screw-capped plastic cryovials. The product is provided frozen on dry ice and should be stored at -20°C or colder immediately upon arrival. Freeze-thaw cycles should be minimized.

Citation:

Acknowledgment for publications should read "The following reagent was obtained through BEI Resources, NIAID, NIH: SARS-Related Coronavirus 2, Wuhan-Hu-1 Spike-Pseudotyped Lentiviral Kit, NR-52948."

Biosafety Level: 1

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. [Biosafety in Microbiological and Biomedical Laboratories](#). 5th ed. Washington, DC: U.S. Government Printing Office, 2009; see www.cdc.gov/biosafety/publications/bmbl5/index.htm.

Disclaimers:

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at www.beiresources.org.

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC® nor the U.S. Government makes any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC® nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. ATCC® and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, ATCC®, their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

Use Restrictions:

This material is distributed for internal research, non-commercial purposes only. This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale.

NR-52948 is claimed in U.S. Patent number 8,008,006 and European Patent number 1341808 and the continuations, continuations-in-part, re-issues and foreign counterparts thereof.

To obtain a license for commercial use, users may contact Fred Hutchinson licensing via email at mta@fredhutch.org.

References:

1. Bloom, J. and A. Balasz, Personal Communication.
2. Crawford, K. H. D., et al. "Protocol and Reagents for Pseudotyping Lentiviral Particles with SARS-CoV-2 Spike Protein for Neutralization Assays." *Viruses* 12 (2020): E513. PubMed: 32384820.
3. Matz, M. V., et al. "Fluorescent Proteins from Nonbioluminescent *Anthozoa* Species." *Nat. Biotechnol.* 17 (1999): 969-973. PubMed: 10504696.
4. Murphy, G. J., et al. "Exogenous Control of Mammalian Gene Expression via Modulation of Translational

Termination." *Nat. Med.* 12 (2006): 1093-1099. PubMed: 16892063.

5. Hulswit, R. J. G., C. A. M. de Haan and B.-J. Bosch. "Coronavirus Spike Protein and Tropism Changes." *Adv. Virus Res.* 96 (2016): 29-57. PubMed: 27712627.

ATCC® is a trademark of the American Type Culture Collection.

